





EcoBino LED

BE211A



EDUCATIONAL / HOBBY MICROSCOPE (BE211A ECO-BINO-LED)

The popular choice of microscope for grade schools, home schooling students and microscopic enthusiasts who love to see things up close and personal. The Eco-Bino compound microscope is designed with the features of minimalist microscopy in mind.

Binocular head at 45°
360° rotating head
WF10x/18mm eyepieces
Revolver with achromatic objectives 4x,10x, 40x and 100x oil
Condenser with Iris diaphragm
LED illumination

Popular starter microscope for the following:

Bee keepers
Grade school teachers
Home schooling students
College/university juniors
Microscopy hobbyists
Hobby botanists